SBURNTEC® FIRST AID FOR BURNS

North American Rescue® is proud to lead the way with new technology as we introduce the **advanced BurnTec**® **first aid dressings** to our customers. BurnTec® is **FDA cleared** and **CE marked**, and it is currently being utilized in European countries. This hydrogel dressing is different from other burn first aid dressings due to its **stable-sheeting hydrogel structure**. The dressing has been strengthened with a medicalquality non-woven fabric which makes it easy to handle. This stable structure eliminates the risk of spilling or damaging the cooling hydrogel. In other words, the hydrogel stays over the burn area and prevents the hydrogel from **spilling** off body surface area and into body cavities. It also **eliminates the** risk of choking during oxygen therapy.

The constitution of the BurnTec® dressing incorporates an **aqueous composition of natural and synthetic polymers** cross-linked by a beam of electrons.

The BurnTec® cools the injured area and ensures a soothing effect for up to 24 hours. The sterile nonstick dressing is absorbent, trapping wound secretion and bacteria within the gel. The dressing is permeable to water vapor and oxygen but impermeable to bacteria which helps protect the wound from external contamination while

allowing it to breathe. The material is biocompatible and does not cause allergic reactions.

BurnTec® is a modern hydrogel dressing used for a variety of skin injuries ranging from first, second, and third-degree burns (thermal, chemical, & electric). This dressing can also be utilized to provide relief from other skin injuries such as sunburns, abrasions, bruises, swelling, or

Immediate
application of
the dressing
reduces
the risk
of hypertrophic
scars. The

insect bites.

properties allow

non-stick

for painless replacement of the dressing and does not destroy granulation or new epithelium.

BurnTec® dressings are now included

in the **U.S. Marine Corps IFAK**

& the **U.S. Navy Authorized Medical Allowance List** and have received **international recognition**for breakthrough innovation at the **42**nd **World Exhibition of Inventions and Innovation Brussels Eureka**, **International Exhibition of Inventions**(IENA) in Nurnberg and the **International**

Exhibition Genius in Budapest.

For the treatment of:


1st, 2nd & 3rd degree Burns
THERMAL CHEMICAL ELECTRIC

Sunburns

Abrasions & skin irritations

Bruises & SwellingINSECT BITES

Other skin injuries



Thanks to thermal conduction parameters as well as high thermal capacity of hydrogel, heat from the skin surface is transferred by conductive transfer.





NORTH AMERICAN RESCUE

www.NARescue.com • 888.689.6277



NOW in the U.S. Marine Corps IFAK

& the **U.S. Navy Authorized Medical Allowance List**

Breathable, moist wound environment that promotes healing & protects from bacteria

Easy-to-handle, stable **non-stick** Hydrogel leaves no residue behind

Provides cool, soothing relief up to 24 hours for a wide range of skin injuries



BurnTec® Burn Dressing Kit Convenient array of BurnTec® dressings in a convenient carrying case. Kit includes 1 each of the following sizes:

2 in. x 2 in. 2.5 in. x 5 in. 4 in. x 4 in. 5 in. x 5 in. 10 in. x 10 in. (Face) 8 in. x 16 in. 12 in. x 16 in. (Face)

Kit Dimensions: L 12 in. x W 15.5 in. x D 3.5 in. Wt: 6 lb 4 oz

Burn Tec® Burn Dressing Kit Item# 80-0361NSN# 6510-01-630-9177